


3 Parts Infrared Sauna Blanket

INSTRUCTION



Dear customer,

Thanks for your choosing our products.

For the best capability and safety, please read the manual before operation.

Please reserve the manual for your reference.

CONTENTS

Specialty	1
Control panle	2
Operation Procedures	3
Technology Parameter.....	4
Packing list.....	5
Important Safety Notice.....	6




Specialty:

The unique, separable, three temperature zone structure makes it possible for customized program for reducing fat at abdomen, waist, and leg. Other features include safe voltage independent control for three temperature zones, safeguarding device, and reinforced waterproof etc. the carpet asks for very little floor space for installation, and is easy to carry with or operate. It is especially suitable for beauty salon and beauty shop that has small floorage. The product can adjust the temperature at the range of 40°C-85°C and the blanket voltage is 36v.

Caution of interrelated products

- 1.The products that you use together with the machine should be with the high quantity and safe function.
- 2.Don't use the products with the element of alcohol, acid and alkali.
- 3.It's necessary to use the moist and collagen products to protect the skin.
- 4.After the treatment, please use the sun block cream every day.It will be helpful if use the essence gel, but please use the gel with anti- irritability and diminish inflammation function

Control Panel



1. Power
2. Working Time Adjust
3. Working Time Indicator

Output 1

4. Start/pause
5. Temperature Adjust(40~85°C degrees centigrade)
6. Temperature Indicator

Output 2

7. Start/pause
8. Temperature Adjust(40~85°C degrees centigrade)
9. Temperature Indicator

Output 3

10. Start/pause
11. Temperature Adjust(40~85°C degrees centigrade)
12. Temperature Indicator

Operation Procedures:

1. Plug the attachment cords to outlet accordingly;
2. Switch the power after checking (indicator on)
3. Press the 'power' button. The preset temperature is 50°C and instant temperature is displayed after 10 seconds.
4. Adjuster the timer and Set the fat- dissolution temperature by pressing the adjustment button;
Note: in order to rapid melting, the temperature can be higher. But, during the treatment, it should be lower and adjust according to the client's request.
5. Press the Start /Pause key;
6. Power off and disconnect when finish;

ATTENTION:

1. It is helpful to achieve desired effect in assist of effective slimming products as well as point massage
2. Envelope necessary region by fresh preserver before dressed in slimming cloth can maintain cloth clean.

Malfunction management:

Before service, please read the manual carefully, If abnormal conditions continue after check, please contact with the supplier.

You can submit an application to the supplier during the guarantee period. Repair is for free in the period as well as correct use according to the manual.

If exceed guarantee period, you could negotiate with the supplier. If the product can be recovered to normal situation, stipendiary repair is accepted as requested.

Technology Parameter

Voltage: AC110V 220V 230V 240V

Frequency: 50Hz 60Hz

Power: 300W

Packing size(table style): 54×40×23 cm G.W.: 10.8 kg

Packing size(with stand): 103×50×41 cm G.W.: 21.4KG

Before using this machine for the first time, read the operating instruction attentively.

Caution:

Please read the introduction carefully. When malfunction occurs. Please contact with the supplier. Do not open it by yourself to be protected. keep out of water to avoid damage or danger when cleansing the machine.

Packing list

Item Name	Unit	QTY	Mark	
Mainframe	Set	1		
Power wire	Strip	1		
Blanket	Piece	1		
Fuse	Piece	2		
Salver	Piece	1	Attched with stand	
Instruction	Piece	1		

Any change is subjected to manufacturer without notice.

Important Safety Notice

1. Please don't disassemble the equipment or try to do some other operation, which have been indicated in our instruction. All the repairing work should be done by our admitted professional personnel
2. Please don't assemble or operate this equipment while you are near the water or your hands are wet. Notice that don't pour the liquid on it.
3. Please pull out the plug from the socket when you met such conditions as below. And keep in touch with professional maintenance Station.
 1. The equipment touched liquid
 2. The equipment has abnormal smell, smoking or loud and strong noise.
 3. Cable is broken
 4. The equipment is fallen down and broken
4. Please don't put anything on the cables, don't place the equipment at any place where someone can step on the cable.
5. Please pull out the power cable after using it for your safety
6. Don't put anything into the equipment's jack, otherwise it may cause fire or touch the electricity. If something enters the equipment, please don't pick it out by yourself, contract distributor or our company immediately
7. Please don't put the equipment on the anything that can't stand steadily, avoid it turns over and damage.
8. Our company has the right to modify the information in this introduction before notify! And keep the explanations right for the above.